**SAFETY INSTRUCTIONS FOR  
ИНСТРУКЦИИ ЗА БЕЗОПАСНОСТ ЗА  
BEZPEČNOSTNÍ POKYNY TÝKAJÍCÍ SE  
SIKKERHEDSINSTRUKTIONER TIL  
SICHERHEITSHINWEISE FÜR  
INSTRUCCIONES DE SEGURIDAD PARA  
OHUTUSEESKIRJAD  
ΟΔΗΓΙΕΣ ΓΙΑ ΤΗΝ ΑΣΦΑΛΕΙΑ ΓΙΑ  
INSTRUCTIONS DE SÉCURITÉ POUR  
SIGURNOSNE UPUTE ZA  
ISTRUZIONI DI SICUREZZA PER  
ÖRYGGISLEIÐBEININGAR FYRIR  
DROŠĪBAS NORĀDĪJUMI  
SAUGOS NURODYMAI, SKIRTI  
BIZTONSÁGI UTASÍTÁSOK  
VEILIGHEIDSVOORSCHRIFTEN VOOR  
SIKKERHETSINSTRUKSER FOR  
INSTRUKCJE BEZPIECZEŃSTWA DLA  
INSTRUÇÕES DE SEGURANÇA PARA  
INSTRUCȚIUNI DE SIGURANȚĂ PENTRU  
BEZPEČNOSTNÉ POKYNY PRE  
VARNOSTNA NAVODILA ZA  
TURVALLISUUSOHJEET  
SÄKERHETSINSTRUKTIONER FÖR**

**ShockMaster 500**

**ShockMaster 300**

Български	2 - 3
Česky	4 - 5
Dansk	6 - 7
Deutsch	8 - 9
English	10 - 11
Español	12 - 13
Eesti keel	14 - 15
Ελληνικά	16 - 17
Français	18 - 19
Hrvatski	20 - 21
Italiano	22 - 23
Íslenska	24 - 25
Latviešu	26 - 27
Lietuvių kalba	28 - 29
Magyar	30 - 31
Nederlands	32 - 33
Norsk	34 - 35
Polski	36 - 37
Português	38 - 39
Română	40 - 41
Slovenščina	42 - 43
Slovensky	44 - 45
Suomi	46 - 47
Svenska	48 - 49

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## 1. Installation

### 1.1 On receipt of the appliance

Check if the appliance and the accessories were not damaged during the shipment. In the event of any damage and/or malfunction notify your supplier immediately (= within 3 working days after the delivery) by means of telephone, fax, e-mail or letter. In the case of damage, the unit may not be operated.

### 1.2 Connecting and starting

Carefully follow the instructions step by step for the installation of the ShockMaster 500.

The ShockMaster must be installed and commissioned in accordance with the EMC guidelines detailed in the EMC guidance and manufacturer's declaration section.

Place the appliance horizontally and stable.

Check the supporting surface has sufficient carrying capacity to avoid equipment damage.

The device must be positioned in such a way that the plug at one end of the detachable power cord can be easily removed from the outlet socket respectively the mains entry socket.

Do not put the appliance in direct sunlight and do not put it above a heat source.

No objects may be put on the appliance.

Make sure the ventilation holes in the housing of the ShockMaster and slots in the trolley are not covered during operation.

It is essential to prevent fluids entering the unit or its hoses.

### 1.3 Functional test

During production the appliance is thoroughly checked for electrical safety.

Check if the display is functioning properly. If this is not the case, do not continue to use the device and inform the repair services of your local GymnaUniphy dealer.

### 1.4 Transport of the appliance

The way it is packed in its original transport packaging, the appliance is suitable for transportation and storage within the following limits:

- Temperature from +0 °C to +60 °C;
- Relative humidity from 5 % - 95 %, inclusive condensation;
- atmospheric pressure from 50 kPa - 106 kPa.

If it is necessary to dispatch the appliance, we recommend that you use the original packaging to do so.

### 1.5 Disposal

Disposal of the device and its components must be carried out in accordance with national waste disposal regulations.

Used cleaning / disinfecting agents must be disposed of according to the user instructions of these agents.

### 1.6 When no printed manual is provided

Upon request and without additional cost, one paper copy of the manual in one of the available languages will be provided promptly to the user for each purchased device.

## 2. Safety aspects

### 2.1 Electrical safety

To avoid the risk of electric shocks the device must only be connected to an outlet socket with a protective earth terminal with the original power cord provided by the manufacturer.

Do not open the device. Risk of electric shocks!

Verify the ratings on the type label of the compressor match with the supply mains properties before connecting the ShockMaster 500 to the power mains.

Disconnect the unit from the mains before cleaning, disinfection or maintenance work.

## **2.2 Explosion prevention**

Do not use the ShockMaster in damp rooms or potentially explosive environments, i.e. in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.

## **2.3 Operational safety**

Nearby portable and mobile HF communications equipment such as cell phones may cause interferences.

Before you start using the ShockMaster for the first time, please make sure you have read and understood all information provided in the user manual and accessory manuals for the detachable components. Thorough knowledge of the information in these manuals will enable you to react promptly and effectively in case of malfunctions and failures.

All status and error messages that appear during treatment must always be acted upon immediately.

The ShockMaster is exclusively intended for use by medical specialists and must only be used by suitably qualified and trained medical personnel.

The user is responsible for correctly positioning the hand piece of the ShockMaster and determining where the treatment zone of the patient is located.

Prior to every new treatment the user must ask the patient for possible contra-indications and allergy to any ingredient of the coupling gel. Furthermore, the user must educate the patient regarding the actions, side effects and possible residual hazards related to the treatment.

To avoid safety hazards, the device must not be used for applications other than these:

- biomechanical therapy
- myofascial trigger points (MTrP)
- disorder of tendon insertions
- activation of muscle and connective tissue
- acupuncture shockwave therapy

The contra-indications listed below are examples. No claims are made regarding the completeness or unlimited validity of the list.

- coagulation disorders (haemophilia)
- use of anticoagulants, especially Marcumar
- thrombosis
- tumour diseases, carcinoma patients
- pregnancy
- polyneuropathy in case of diabetes mellitus
- acute inflammations / pus focus in the target area
- open epiphyseal discs
- cortisone therapy up to 6 weeks before first treatment
- patient with a pacemaker
- prosthesis
- osteoporosis
- infected wounds
- large nerves and vessel cavities that contain air (lung, intestines, ...)
- risk of haemorrhaging
- cardiac region
- open scar
- vertebrae, spinal column or head

Side effects generally disappearing within 5 to 10 days.

- swelling, reddening, haematomas
- petechiae, pain, irritation of the periosteum
- skin lesions after previous cortisone therapy
- cardiac arrhythmias

NOTE:

The operating ambient temperature for the appliance is between +10°C and +40°C.

## **2.4 Medical Device Directive**

The appliance meets the most recent adapted essential requirements stated in the Medical Device Directive issued by the European Commission (93/42/EEC).

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