User Manual

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
Manual do Utilizador

Manuale d'uso

VACO 200

English	3
Nederlands	29
Français	55
Deutsch	81
Español	109
Português	135
Italiano	161

© 2007, GymnaUniphy N.V.

All rights reserved. Nothing from this publication may be copied, stored in an automated data file, or made public, in any form or in any way, be it electronically, mechanically, by photocopying, recordings or in any other way, without prior written permission from GymnaUniphy N.V.

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User Manual Vaco 200

Device for the application of electrotherapy using vacuum electrodes

Manufacturer	GymnaUniphy N.V.
Main office	Pasweg 6A B-3740 BILZEN
Telephone	+ (32) (0)89-510.510
Fax	+ (32) (0)89-510.511
E-mail	info@gymna-uniphy.com
Website	www.gymna-uniphy.com

Version 1.0
June 2007

Abbreviations

EMC	Electromagnetic Compatibility
EL	Electrode
ET	Electrotherapy

Symbols on the equipment

Read the manual.

Symbols in the manual

Warning or important information.

TABLE OF CONTENTS

1	SAFETY	7
1.1	PURPOSE	7
1.2	SAFETY INSTRUCTIONS	7
1.3	MEDICAL DEVICES DIRECTIVE	9
1.4	LIABILITY	9
2	INSTALLATION	10
2.1	RECEIPT	10
2.2	PLACING AND CONNECTION	10
2.3	USE IN COMBINATION WITH ANOTHER DEVICE	10
2.4	TRANSPORT AND STORAGE	10
2.5	RESELLING	10
3	DESCRIPTION OF THE EQUIPMENT	11
3.1	VACO 200 AND STANDARD ACCESSORIES	11
3.2	COMPONENTS OF VACO 200	12
4	OPERATION	13
4.1	PERFORMING ELECTROTHERAPY WITH THE VACO 200	13
4.2	CONNECTING THE ELECTRODES	14
4.3	THE WATER RESERVOIR IS FULL	16
4.4	REFERENCE	16
5	INSPECTIONS AND MAINTENANCE	17
5.1	INSPECTIONS	17
5.2	MAINTENANCE	18
6	MALFUNCTIONS, SERVICE AND GUARANTEE	20
6.1	MALFUNCTIONS	20
6.2	SERVICE	21
6.3	GUARANTEE	21
6.4	TECHNICAL LIFE TIME	22
7	TECHNICAL INFORMATION	23
7.1	GENERAL	23
7.2	ENVIRONMENTAL CONDITIONS	23
7.3	TRANSPORT AND STORAGE	23
7.4	STANDARD ACCESSORIES	24
7.5	OPTIONAL ACCESSORIES	24

8	APPENDICES	25
	8.1 TECHNICAL SAFETY INSPECTION	25
	8.2 EMC DIRECTIVE	26
	8.3 DISPOSAL	26
	8.4 LITERATURE	26
9	INDEX	27
10	REFERENCE EMC DIRECTIVE	189
	10.1 GUIDANCE AND DECLARATION	189

1 SAFETY

1.1 Purpose

The Vaco 200 is intended solely for medical applications. The Vaco 200 is a vacuum suction device for fast and efficient electrode placement. The device is suited for continuous use.

1.2 Safety instructions

1.2.1 General

- Only qualified people who are trained in the application of the therapies may use the appliance.
- Only a technician authorised by GymnaUniphy N.V. may open the equipment or the accessories.
- Follow the instructions and directions in these user instructions.
- Place the equipment on a horizontal and stable base.
- Do not place the equipment in the sun or above a heat source.
- Do not use the equipment in a damp area.
- Do not let any liquid flow into the equipment.
- Do not disinfect or sterilise the equipment. Clean the equipment with a dry or moistened cloth. *See §5.2.*
- The 'Directive on Medical Devices' from the European Commission (93/42/EEG) requires that safe devices are used. It is recommended to perform a yearly technical safety inspection. *See §5.1.1.*
- Disconnect the equipment from the patient before switching the power on or off.
- Be careful, in combination with an electrotherapy device the output of the Vaco 200 can have physiological effects. See the 'Safety instructions' in the user manual of the electrotherapy device.
- Always keep these user instructions with the equipment.

1.2.2 **Electrical safety**

- Only use the equipment in an area with facilities that meet the applicable legal regulations.
- Connect the equipment to an outlet with a protective earth terminal. The outlet must meet the locally applicable requirements for medical areas.

1.2.3 **Prevention of explosion**

- Do not use the equipment in an area where combustible gases or vapours are present.
- Switch off the equipment when it is not used.

1.2.4 **Electro Magnetic Compatibility**

- Medical electrical equipment requires special precautions for Electro Magnetic Compatibility (EMC). Follow the instructions for the installation of the equipment. *See §2.*
- Do not use mobile telephones or other radio, shortwave, or microwave equipment in the vicinity of the equipment. This kind of equipment can cause disturbances.
- When the device operates in the vicinity of equipment other than those mentioned in §2.3, verify that the Vaco 200 operates normally.
- Only use the accompanying accessories that are supplied by GymnaUniphy. *See §7.4 and §7.5.* Other accessories can lead to an increased emission or a reduced immunity.

1.2.5 **Accessories**

- Check the electrode cables and the electrodes at least once a month. Check whether the insulation is still intact. *See §5.1.*
- Use always demineralized water with vacuum electrodes to avoid lime deposits in the watertank, the tubes and the sponges.
- Use only moisten sponges. Too dry sponges can cause burning of the skin under the electrodes.
- Do not use vacuum electrodes with DC currents. The DC currents cause damage to the vacuum cups by ionization.

1.3 Medical Devices Directive

The device complies with the essential requirements of the Medical Device Directive of the European Committee (93/42/EEC) as most recently changed.

1.4 Liability

The manufacturer cannot be held liable for injury to the therapist, the patient or third parties, or for damage to or by the equipment used, if for example:

- an incorrect diagnosis is made.
- the equipment or the accessories are used incorrectly.
- the user instructions are wrongly interpreted or ignored.
- the equipment is badly maintained.
- maintenance or repairs are performed by people or organisations that are not authorised to do so by GymnaUniphy.

Neither the manufacturer nor the local GymnaUniphy dealer can be held liable, in any way whatsoever, for the transfer of infections via accessories.

2 INSTALLATION

2.1 Receipt

1. Check whether the equipment has been damaged during transport.
2. Check whether the accessories are intact and complete. *See §7.4.*
 - Inform your supplier of any damage or defects by no later than within 3 working days after receipt. Report the damage by telephone, fax, e-mail or letter.
 - Do not use the equipment if it is damaged or defective.

2.2 Placing and connection

1. Place the equipment on a horizontal and stable base.
 - Do not place the equipment in the sun or above a heat source.
 - Do not use the equipment in a wet area.
2. Check whether the voltage that is stated on the rear of the equipment corresponds with the voltage of your mains supply. The equipment is suited for a nominal mains voltage from 100 V to 240 VAC / 50-60 Hz.
3. Connect the device to an outlet with protective earth terminal.

2.3 Use in combination with another device

The Vaco 200 can be used in combination with:

- The Combi200.
- The Combin200L.
- The Duo200.

2.4 Transport and storage

Take account of the following matters if the equipment has to be transported or stored:


- Transport or store the equipment in the original packaging.
- The maximum period for transport or storage is: 15 weeks.
- Temperature: -20 °C to +60 °C.
- Relative humidity: 10% to 100%.
- Atmospheric pressure: 200 hPa to 1060 hPa.

2.5 Reselling

This medical equipment must be traceable. The equipment and some accessories have a unique serial number. Provide the dealer with the name and address of the new owner.


3 DESCRIPTION OF THE EQUIPMENT

3.1 Vaco 200 and standard accessories

1. Vaco 200. See §3.2.
2. Power cord
3. Vacuum hose (4 pieces)
4. Connection cable: ET device - Vaco 200 (2 pieces)
5. Sponge for vacuum electrode (4 pieces)
6. Vacuum electrode (4 pieces)

3.2 Components of Vaco 200

- | | |
|---|---|
| 1. Vacuum intensity regulator | for two-ply electrode cable |
| 2. Pulse regulator | 12. LED: vacuum output channel B active |
| 3. LED water reservoir | 13. Output connectors \llcorner_B for vacuum electrodes |
| 4. Indication: Read manual | 14. Push button channel B |
| 5. Indication: Floating patient circuit | 15. Input connector channel B for connection cable |
| 6. LED: vacuum output channel A active | 16. Input connector channel A for connection cable |
| 7. Output connectors \llcorner_A for vacuum electrodes | 17. On/off switch |
| 8. LED: traditional output channel A active | 18. Fuse holder |
| 9. Output connector \llcorner_A for two-ply electrode cable | 19. Connection to mains supply |
| 10. LED traditional output channel B active | 20. Drain hose for water reservoir |
| 11. Output connector \llcorner_B | 21. Air exhaust |
| | 22. Type plate |

4 OPERATION

4.1 Performing electrotherapy with the Vaco 200


1. Connect the Vaco 200 and the electrotherapy device to the mains supply.
2. Connect the input connector A of the Vaco 200 to the output connector A of the electrotherapy device with the connection cable. If necessary, do the same for channel B.
3. Turn the vacuum intensity regulator to the 'Y' position.
4. Use the switch on the rear panel to set the Vaco 200 and the electrotherapy device to on. The LEDs from the traditional channels of the Vaco 200 light up. The LEDs show which output channel is connected with the electrotherapy device.
5. Select the desired electrotherapy program.
6. Connect, adjust and place the electrodes. *See §4.2.1, §4.2.2, §4.2.3 and §4.2.4.*
7. Turn the intensity knob A or B of the ET device to start the electrotherapy and to set the desired intensity.
8. Check the patient's reaction. Repeat this check regularly during the treatment.
9. The electrotherapy device stops the treatment and indicates that the treatment is completed.
10. If vacuum is used, turn the vacuum intensity regulator to 'Y'. For safety, the vacuum stops automatically two minutes after the treatment stops and the vacuum channel LEDs flash. Turn the vacuum intensity regulator to 'Y' to make a new treatment selection possible.
11. Remove the electrodes. If necessary put your finger under the rim of the vacuum electrode for air inled.


4.2 Connecting the electrodes


- Use always demineralized water with vacuum electrodes to avoid lime deposits in the watertank, tubes and sponges. Add a saline solution to improve the electrical conduction.
- Use moist sponges only. Too dry sponges can cause a bad electrical contact and burn the skin.
- Do not use vacuum electrodes with DC current. The DC current causes damage to the vacuum cups by ionization.


4.2.1 Using four vacuum electrodes

1. If present disconnect the vaginal, anal and rectal stimulation probes.
2. Connect the vacuum electrodes to the vacuum hoses.
3. Connect the four vacuum hoses. Select two cables with the same hose colour for each channel.
 1. Connect the red connectors from the vacuum hoses to the output connectors with the red dot.
 2. Connect the black connectors from the vacuum hoses to the output connectors with the black dot.
4. Moisten the round sponges.
5. Put the sponges in the vacuum electrodes.
6. Turn the vacuum intensity regulator to  and regulate the desired suction force.
7. Place the vacuum electrodes on the part of the body that must be treated. The vacuum electrodes stay in place by the underpressure. Too high suction cause patient discomfort.
8. Turn the pulse vacuum suction power regulator to  for patient comfort.
9. Select the stimaton signal on the vacuum output  with the push button. The LED vacuum channel B is lit.


4.2.2 Using two vacuum electrodes

1. If present disconnect the vaginal, anal and rectal stimulation probes.
2. Connect two vacuum electrodes to the vacuum hoses. Select two cables with the same hose colour.
3. Connect the red connector from the vacuum hose to the channel A output connector with the red dot.
4. Connect the black connector from the vacuum hose to the channel A output connector with the black dot.
5. Connect a vacuum hose for a closed vacuum system to make a short-circuit across the channel B output connectors.
6. Moisten the round sponges.
7. Put the sponges in the vacuum electrodes.
8. Turn the vacuum intensity regulator to  and regulate the desired suction force.
9. Place the vacuum electrodes on the part of the body that must be treated. The vacuum electrodes stay in place by the underpressure. Too high suction cause patient discomfort.
10. Turn the pulse vacuum suction power regulator to  for the patient comfort.

4.2.3 Using traditional electrodes


1. Prepare the traditional electrodes. *See the User Manual of the ET device.*
2. Connect the two-ply cable to the corresponding output connector ∇_A or ∇_B of the Vaco 200.

4.2.4 Using a combination of vacuum and traditional electrodes

Use channel A for the stimulated vacuum electrodes and use channel B for the traditional electrodes.

1. Connect two vacuum electrodes to channel A and short-circuit channel B. See §4.2.2.
2. Prepare the traditional electrodes. See the *User Manual of the ET device*.
3. Connect the two-ply cable to the output connector Ψ_B of the Vaco 200.

4.3 The water reservoir is full

1. The LED of the water reservoir  light up. Finish the present treatment. After you turn the vacuum off, you can not restart the vacuum.
2. Empty the water reservoir with the drain hose. The start of a vacuum treatment is possible.
3. Clean the water reservoir. See §5.2.4.

4.4 Reference

4.4.1 Indications

The Vaco 200 can be used for the pathologies from the indication list of the electrotherapy device. See the *User Manual of the ET device*.

4.4.2 Contra indications

- The contra indications from the electrotherapy device. See the *User Manual of the ET device*.
- Internal infections
- Hemorrhagic risk in the part of the body where the electrodes are to be placed.

5 INSPECTIONS AND MAINTENANCE

5.1 Inspections

Component	Check	Frequency
Electrode cables and electrodes	Damage Insulation intact	At least 1x per month
Equipment	Technical safety inspection. See §5.1.1.	At least 1x per year

5.1.1 Technical safety inspection

The 'Directive on Medical Devices' from the European Commission (93/42/EEG) requires that safe devices are used. It is recommended to perform a yearly technical safety inspection. If the legislation in your country or your insurer prescribes a shorter period, you must adhere to this shorter period.

- Only a technician authorised by GymnaUniphy N.V. may open the equipment or the accessories.
- The inspection may only be performed by a suitably qualified person. In some countries this means that the person must be accredited.

Inspection points

The technical safety inspection include the following tests:

1. Test 1: General: Visual inspection and check on the operating functions
2. Test 2: Electrical safety inspection: measurement of the earth leakage current and patient leakage current according to DIN/VDE 0751-1 ed. 2.0.

Inspection result

1. A registration must be maintained of the technical safety inspections. Use the inspection report in the appendix for this purpose. See §8.1.
2. Copy this appendix.
3. Complete the copied appendix.
4. Keep the inspection reports for at least 10 years.

The inspection is successful if all inspection items are passed.
Repair all faults on the equipment before the equipment is put back into operation.

By comparing the registered measurement values with previous measurements, a possible slowly-deteriorating deviation can be ascertained.

5.2 Maintenance

Component	Check	Frequency
Vaco 200 unit	Cleaning. See §5.2.1.	If necessary
Vacuum electrodes	Cleaning. See §5.2.2.	After every treatment
Sponges	Cleaning. See §5.2.3.	After every treatment
Vacuum hoses and water reservoir	Cleaning. See §5.2.4.	Weekly

Accessories that come in contact with the body of the patient must be washed with pure water after the disinfection to prevent allergic reactions.

5.2.1 Cleaning the Vaco 200 unit

1. Switch the equipment off.
2. Clean the device with a lightly moistened soft cloth.

No liquid may enter the device.

5.2.2 Cleaning the vacuum electrodes

1. Clean the vacuum electrodes (metal electrodes and rubber suction cups) in a non-aggressive soap solution or in a 70% alcohol solution.
2. Rinse the vacuum electrodes thoroughly with water.
3. Dry the vacuum electrodes.

5.2.3 Cleaning the sponges for vacuum electrodes

1. Rinse the sponges thoroughly with water or clean the sponges with a 70% alcohol solution.
2. Rinse the sponges thoroughly with water.

5.2.4 *Cleaning the vacuum hoses and the water reservoir*

1. Empty the water reservoir with the drain hose.
2. Connect the vacuum hoses.
3. Put the ends of the vacuum hoses in a 70% alcohol solution.
4. Turn the vacuum unit on.
5. Suck the liquid up until the LED of the water reservoir lits.
6. Turn the vacuum unit off.
7. Empty the water reservoir.
8. Repeat steps 3 until 7 with pure water.

6 MALFUNCTIONS, SERVICE AND GUARANTEE

6.1 Malfunctions

Component	Problem	Solution
Vaco 200 unit	Equipment cannot be switched on	See §6.1.1.
	Equipment does not react to commands	See §6.1.3.
Vacuum electrodes	Contamination by ionization	See §6.1.4.
Sponges	Furring	Replace the sponges
	Bad conduction	Replace the sponges

6.1.1 *Equipment cannot be switched on*

1. Check if the mains voltage has failed.
2. Check if the main switch is switched on ("I").
3. Check if the power cord and the fuses are in order. If necessary, replace the fuses. See §6.1.2.
4. Contact your dealer if the equipment still cannot be switched on.

6.1.2 *Replacing the fuses*

1. Switch the main switch off ("O").
2. Unplug the power cord from the equipment.
3. Pull the fuse holder carefully out of the equipment. If necessary, use a screwdriver.
4. Replace both fuses. If necessary, order new fuses from your dealer.
5. Install the fuse holder and plug in the power cord.
6. Switch the main switch on again ("I").

6.1.3 *Equipment does not react to commands*

The safety system of the equipment has ascertained a fault. You cannot continue to work.

1. Disconnect the connection to the patient.
2. Switch the main switch off ("O").
3. Wait 5 seconds and switch the main switch on again ("I").
4. Contact your dealer if the fault reappears.

6.1.4 Remove the contamination from the vacuum electrodes

1. Clean the vacuum electrodes. See §5.2.2.
2. Use steel wool or sandpaper for metal with fine grains ('P 400' or higher) to remove the contamination.
3. Replace the vacuum electrodes if the contamination is still present.

6.2 Service

Only a technician authorised by GymnaUniphy N.V. may open the equipment or the accessories to perform repairs. The equipment does not contain any components that may be replaced by the user.

Service and guarantee are provided by your local GymnaUniphy dealer. The conditions of delivery of your local GymnaUniphy dealer apply. If you have qualified technical personnel that are authorised by GymnaUniphy to perform repairs, your dealer can provide diagrams, spare parts lists, calibration instructions, spare parts and other information on request, for a fee.

6.3 Guarantee

GymnaUniphy and the local GymnaUniphy dealer declares itself to be solely responsible for the correct operation when:

- all repairs, modifications, extensions or adjustments are performed by authorised people.
- the electrical installation of the relevant area meets the applicable legal regulations.
- the equipment is only used by suitably qualified people, according to these user instructions.
- the equipment is used for the purpose for which it is designed.
- maintenance of the device is regularly performed in the way prescribed. See §5.
- the technical life time of the equipment and the accessories is not exceeded.
- the legal regulations with regard to the use of the equipment have been observed.

The guarantee period for the equipment is 2 (two) years, beginning on the date of purchase. The date on the purchase invoice acts as proof. This guarantee covers all material and production faults. Consumables, such as sponges, are excluded.

This guarantee does not apply to the repair of defects that are caused:

- by incorrect use of the equipment,
- by an incorrect interpretation or not accurately following these user instructions,
- by lime deposits,
- by carelessness or misuse,
- as a consequence of maintenance or repairs performed by people or organisations that are not authorised to do so by the manufacturer.

6.4 Technical life time

The expected life time of the equipment is 10 years, calculated from the date of manufacture. See the type plate for this information.

In so far as possible, GymnaUniphy will supply service, spare parts and accessories for a period of 10 years from the date of manufacture.

7 TECHNICAL INFORMATION

7.1 General

Dimensions Vaco 200 (w x h x d)	267 x 95 x 270 mm
Weight Vaco 200	3,5 kg
Weight including accessories	4,6 kg
Mains voltage	100 - 240 VAC, 50 - 60 Hz
Maximum power consumption	30 VA
Safety class	Class I (earthed socket required)
Insulation	Type BF (floating patient circuit)
Fuses	2 x T2AL250V
Volume water reservoir	± 180 ml
Working pressure continuous vacuum	38 - 320 hPa
Working pressure pulsation vacuum	58 - 480 hPa
Vacuum rhythm	1,5/1,5 - 1,5/4,5 s (on/off time)

7.2 Environmental conditions

Temperature	+10 °C to +40 °C
Relative humidity	30% to 75%
Atmospheric pressure	700 hPa to 1060 hPa

7.3 Transport and storage

Transport weight	5,6 kg
Storage temperature	-20 °C to +60 °C
Relative humidity	10% to 100%, including condensation
Atmospheric pressure	200 hPa to 1060 hPa
Transport classification	Single piece by mail

The transport and storage specifications apply to equipment in the original packaging.

7.4 Standard accessories

	Quantity	Description	Art. no.
	1	Power cord ¹	100.689
	2	Connection cable: ET device - Vaco 200	102.032
	1	Vacuum hose dark grey (per 2 pcs: black/red connector)	102.801
	1	Vacuum hose light grey (per 2 pcs: black/red connector)	102.800
	2	Vacuum electrode - 60 mm (per 2 pcs)	114.688
	1	Sponge for vacuum electrode - 60 mm (per 4 pcs)	114.689
	1	User manual Vaco 200 EN/NL/FR/DE/ES/IT	311.025
	1	Safety instructions	323.011

¹ This power cord has a CEE 7/7 type plug. For countries with other outlets, a different power cord with the appropriate plug is supplied.

7.5 Optional accessories

	Quantity	Description	Art. no.
	1	Clips for 2 vacuum tubes	112.457
	2	Vacuum electrode - 90mm	114.686
	4	Sponge for vacuum electrode - 90mm	114.687

Article numbers can change in the course of time. Check the article numbers in the most recent catalogue or ask your dealer.
The drawings are merely indicative, no rights can be derived from them.

8 APPENDICES

8.1 Technical safety inspection

Vaco 200 with serial number is / is not¹ in good working order		
Location:	Inspection performed by: Name:	Owner: Name:
Date:	Initials:	Initials:

¹ Cross out what does not apply.

If a specific test does not apply to this equipment, place a mark in the NA (not applicable) column.

8.1.1 Test 1: General

	Yes	No	NA
1. The results of earlier safety inspections are available.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
2. The logbook is present.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
3. The type plate and the supplier's label are legible.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
4. The housing, adjusting knobs and keys are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
5. The power connection and power cord are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6. The output connectors are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
7. The electrode connectors and cables are undamaged.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8. All the operation functions work.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
9. The measured output signals are correct.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
10. All the alarm functions work.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8.1.2 Test 2: Electrical safety test (VDE 0751)

	Yes	No
1. The resistance of the safety earth is less than 0.2 Ω	<input type="checkbox"/>	<input type="checkbox"/>
2. The housing leakage current is less than 1000 μA	<input type="checkbox"/>	<input type="checkbox"/>
3. The patient leakage current is less than 5000 μA	<input type="checkbox"/>	<input type="checkbox"/>

Notes:

8.2 EMC directive

The tabel at the last pages of the user manual shows the information about the EMC properties of the equipment. See §10.1. Because this information is intended for technicians, the information is given in English.

8.3 Disposal

Take account of the following environmental aspects when disposing of the equipment and the accessories:

- The basic device, the cables and the electrodes fall under small chemical waste (or electronic waste). These components contain lead, tin, copper, iron, various other metals and various plastics, etc. Consult the applicable national regulations.
- Sponges contain only organic material and do not require any special processing.
- Packaging materials and manuals can be recycled. Deliver them to the appropriate collection points or include them with the normal household waste. This depends on the local organisation of the waste processing.

8.4 Literature

A literature list can be sent on request. Please contact GymnaUniphy.

INDEX**A**

Abbreviations 4
 Accessories 24
 safety 8

C

Cleaning 18
 Connection 10
 Contamination
 remove 21
 Contra indications 16

D

Device
 cleaning 18
 Disposal 26

E

Electrical safety 8
 Electrodes
 connection 14
 Electrotherapy
 perform 13
 EMC 8
 EMC directive 26
 Environmental conditions 23

G

Guarantee 21

I

Indications 16
 Inspections 17
 Installation 10

L

Liability 9

M

Maintenance 18
 Malfunctions 20
 Medical Devices Directive 9

P

Placing 10
 Prevention of explosion 8
 Purpose 7

R

Replacing the fuses 20
 Reselling 10

S

Safety 7
 instructions 7
 technical inspection 17, 25
 Service 21
 Sponges
 cleaning 18
 Storage 10
 conditions 23

T

Technical information 23
 Technical life time 22
 Transport 10
 conditions 23

V

Vacuum electrode
 cleaning 18
 Vacuum electrodes
 place 14
 Vacuum hoses
 cleaning 19

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10 REFERENCE EMC DIRECTIVE

10.1 Guidance and declaration

Guidance and manufacturer's declaration - electromagnetic emissions		
The 200-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a 200-series device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The 200-series devices use RF energy only for their internal function. Therefore, their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
	Class B	
Harmonic emissions IEC 61000-3-3	Class B	The 200-series devices are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The 200-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a 200-series device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact / ± 8 kV air No loss of performance	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV power / ± 1 kV I/O No loss of performance	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV diff. / ± 2 kV comm. No loss of performance	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 0,5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 sec	$U_T - 100\%$ (0,5 period) No loss of performance $U_T - 60\%$ (5 periods) No loss of performance $U_T - 30\%$ (25 periods) No loss of performance $U_T - 100\%$ (5 seconds) Device resets to a safe state. (60601-1 § 49.2)	Mains power quality should be that of a typical commercial or hospital environment. If the user of a 200-series device requires continued operation during power mains interruptions, it is recommended that the 200-series device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not applicable	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level			

Guidance and manufacturer's declaration - electromagnetic immunity			
The 200-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a 200-series device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} AM 1 kHz 80% 150 kHz to 80 MHz	10 V.....0,15-80 Mhz 26 V.....6,78 Mhz 32 V.....13,56 Mhz 14 V.....27,12 Mhz 27 V.....40,68 Mhz	Portable and mobile RF communications equipment should be used no closer to any part of a 200-series device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0,35\sqrt{p}$ $d = 0,46\sqrt{p}$ $d = 0,38\sqrt{p}$ $d = 0,86\sqrt{p}$ $d = 0,44\sqrt{p}$
Radiated RF IEC 61000-4-3	3 V/m AM 1 kHz 80% 80 MHz to 2,5 GHz	10 V/m..0,08-1,0 Ghz 26 V/m...1,4-2,0 Ghz 30 V/m...433,92 Mhz 30 V/m.....915 Mhz	$d = 0,35\sqrt{p}$ 80 MHz to 800 MHz $d = 0,70\sqrt{p}$ 800 MHz to 2,5 GHz $d = 0,12\sqrt{p}$ $d = 0,23\sqrt{p}$
Radiated RF ENV 50204	3 V/m CW 200 Hz d.c. 50% 895 MHz to 905 MHz	30 V/m.895-905 Mhz	$d = 0,23\sqrt{p}$

Guidance and manufacturer's declaration - electromagnetic immunity

The 200-series devices are intended for use in the electromagnetic environment specified below. The customer or the user of a 200-series device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			<p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 The guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey can be considered. If the measured field strength in the location in which a 200-series device is used exceeds the applicable RF compliance level above, the 200-series devices should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the 200-series device.

b Over the frequency range 150 kHz to 80 MHz, field strengths must be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the 200-series device

The 200-series device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of a 200-series device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 200-series devices as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0,35\sqrt{P}$	80 MHz to 800 MHz $d = 0,35\sqrt{P}$	800 MHz to 2,5 GHz $d = 0,70\sqrt{P}$
0,01	0,04	0,04	0,07

Recommended separation distances between portable and mobile RF communications equipment and the 200-series device

The 200-series device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of a 200-series device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 200-series devices as recommended below, according to the maximum output power of the communications equipment.

0,1	0,11	0,11	0,22
1	0,35	0,35	0,70
10	1,11	1,11	2,21
100	3,50	3,50	7,00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.


Pasweg 6A
B-3740 Bilzen

Tel.: (+32) (0) 89/510.510
Fax: (+32) (0) 89/510.511

www.gymna-uniphy.com
E-mail: info@gymna-uniphy.com

Your dealer / Uw dealer / Votre
distributeur / Ihr Fachhändler / Su
distribuidor / O seu fornecedor / Il
distributore

CE 0344