
EN-Dynamic

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

Copyright:


Enraf-Nonius B.V.
P.O. Box
2600 AV ROTTERDAM
The Netherlands
Tel: +31 (0)15 – 26 98 400
Fax: +31 (0)15 – 25 61 686
info@enraf-nonius.nl
www.enraf-nonius.com

Part number: 1412791_45
3rd edition, Mar 2007

EN-Dynamic

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

NL

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1	Foreword	17
2	Product liability	17
3	Safety instructions	17
3.1	Legal aspects	17
3.2	Warnings	18
3.3	Caution	18
3.4	Note	18
4	Installation	18
5	Intended use	19
6	Contra-indications	19
6.1	Absolute contra-indications	19
6.2	Relative contra-indications	19
7	Operation	19
7.1	EN-Dynamic user interface (with track)	19
7.2	General hardware operation	20
7.3	Seated Leg Press	20
7.4	Squat	21
7.5	Leg Extension	21
7.6	Seated Leg Curl	21
7.7	Abduction	22
7.8	Adduction	22
7.9	Total Hip	22
7.10	Pull Down	23
7.11	Shoulder Press	23
7.12	Chest Press	23
7.13	Rowing	23
7.14	Triceps Dips	24
7.15	Fly	24
7.16	Abdominal Trainer	24
7.17	Low Back Trainer	25
7.18	Upper Body Rotator	25
8	Cleaning and maintenance	25
8.1	Daily maintenance	25
8.2	Weekly maintenance	25
8.3	Periodical maintenance	26
8.4	Annual maintenance	26
9	Finally	26
10	Technical data	26
	Environment conditions for transport and storage	28
10.1	Power supply	28
10.2	Classification	28
10.3	Installation requirements EN-Dynamic	28
10.4	Installation requirements EN-Track	29
10.5	Guidance and manufacturer's declaration –electromagnetic immunity	29
10.6	Guidance and manufacturer's declaration –electromagnetic immunity	30
10.7	Recommended separation distances between portable and mobile RF communications equipment and the EN-Dynamic	30
11	Ordering data	31

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1 Foreword

We very much appreciate your purchase of the EN-Dynamic. The EN-Dynamic offers the possibility to practice in a very flexible way medical training therapy.

When developing and producing EN-Dynamic, much consideration was given to reliability, safety, user comfort and durability. In order to be able to use your EN-Dynamic safely, it is important that you read the instructions before use.

2 Product liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortages of the product.

To the maximum extent permitted by applicable law, in no event will Enraf-Nonius or its suppliers or resellers be liable for any indirect, special, incidental or consequential damages arising from the use of or inability to use the product, including, without limitation, damages for loss of goodwill, work and productivity, or any and all other commercial damages or losses, even if advised of the possibility thereof, and regardless of the legal or equitable theory (contract, tort or otherwise) upon which the claim is based. In any case, Enraf-Nonius's entire liability under any provision of this agreement shall not exceed in the aggregate the sum of the fees paid for this product and fees for support of the product received by Enraf-Nonius under a separate support agreement (if any), with the exception of death or personal injury caused by the negligence of Enraf-Nonius to the extent applicable law prohibits the limitation of damages in such cases.

Enraf-Nonius can not be held liable for any consequence resulting from incorrect information provided by it's personnel, or errors incorporated in this manual and / or other accompanying documentation (including commercial documentation)

The opposing party (product's user or it's representative) shall disclaim Enraf-Nonius from all claims arising from third parties, whatever nature or whatever relationship to the opposing party

3 Safety instructions

Safe execution of the procedures in this manual requires technical experience and general knowledge of safety precautions as well as general knowledge of the human movement and physiological system.

Warnings, Cautions and Notes have been used throughout the manual to bring special matters to the immediate attention of the reader.

- The Warning concerns danger to safety of the user or technician
- The Caution concerns an act that might damage the equipment
- The Note points out a statement deserving more emphasis than the general text, but which does not deserve a Warning or Caution

The sequence of steps in a procedure is important for the safe handling of the equipment and prevention of damage. Therefore never change the sequence of a procedure.

3.1 Legal aspects

This manual and the information herein are the property of Enraf-Nonius B.V. (Rotterdam, The Netherlands).

Enraf-Nonius B.V. disclaims any responsibility for personal injury and/or damage to equipment caused by:

- Negligence or disregard of a Warning or a Caution
- Deviation of any of the described procedures
- Execution of activities that are not described in this manual
- Ignorance of the general safety precautions for handling tools and the use of electricity and pneumatics

- Usage of the device in another manner as specified in this manual.

3.2 Warnings:

- Do not change the load settings during exercising or training: load changes are only allowed in the starting position
- For safety reasons keep sufficient distance from a unit used for training; one should not be within 1 meter distance from the outer borders of a unit while it's used for training
- If pain is experienced during use: stop the training immediately and consult professional medical help
- Use only the black grips to place your hands on.
- Place your feet during training only on the black footsteps. Do not step on the frame.
- EN-Dynamic needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided elsewhere in the usermanual (see technical specifications).
- Portable and mobile Radio Frequency communications equipment can affect EN-Dynamic.
- EN-Dynamic should not be used adjacent to other equipment. If adjacent to other equipment is necessary, the Endynamic should be observed to verify normal operation in the configuration in which it will be used.
- The EN-Dynamic should be operated or used by one person at the time
- Limbs must be kept clear from moving parts
- During exercising hands must be kept clear from the pneumatic cylinder
- Never re-set resistance during exercising
- Keep spectators away from the system when it is in use (at least 1 meter)
- At all times, keep away unsupervised children from the equipment
- A correct biomechanical positioning of the body is necessary
- When the EN-Dynamic is not used in the correct way, injury may occur. Supervision by a knowledgeable (para) medic is therefore obligatory
- The maximum permissible users weight is 135 kg
- Maintain your devices according to the prescribed guidelines (see section 8).

3.3 Caution

- EN-Dynamic is powered by internal or external low voltage, the pneumatic resistance is fed by compressed air. Use only Enraf-Nonius supplied or prescribed spare parts. A wrong voltage or a wrong pressure can result in damage of the unit, or it may result in increased emissions or decreased immunity of the EN-Dynamic.
- The manufacturer of the compressor advises frequent maintenance and the use of special oil lubricants. To prevent malfunction follow the instructions in the service manual of the compressor.
- Take care to prevent the compression of electrical circuits wires and air tubes between the moving parts and the unit itself.

3.4 Note

- The maximum resistance setting for the training equipment will be dependant on the pressure of the compressor.
- In case of any malfunction of your EN-Dynamic unit contact your Enraf-Nonius supplier or Enraf-Nonius B.V. (Rotterdam, The Netherlands). Have repairs, maintenance and installations carried out by an authorized Enraf-Nonius service technician.
- Repairs or maintenance performed by a non Enraf-Nonius authorized service technician and/or the use of (spare) parts that are not originally Enraf-Nonius made, results in a direct loss of any guarantee
-

4 Installation

1. Take the EN-Dynamic and any extra parts out of the packaging and check everything for possible transport damage.
2. Place your device on a flat and hard surface. With an uneven or soft floor you must take measures so that the device cannot rock during use.
3. Read this manual carefully before using the equipment.
4. Mains supply connections must comply with the national requirements regarding medical rooms.
5. Prior to connection of this apparatus to the mains supply, check that the voltage and frequency stated on the type plate correspond with the available mains supply.
6. Avoid exposure to direct sunlight, excessive dust, moisture, mechanical vibrations and shocks.
7. Connection of an external adapter other than the prescribed power supply unit can adversely affect the safety of the user and the functioning of the unit and is therefore not permitted (see technical specifications).

8. The EN-Dynamic must be installed in such a way that its use is not limited by other objects. Keep at least 1 m free around the unit.
9. Do not use the EN-Dynamic in "wet rooms".
- 10.

5 *Intended use*

- Seated Leg Press: Exercising knee and hip muscles
- Squat: Exercising knee and hip muscles
- Leg Extension: Exercising knee extensors
- Seated Leg Curl: Exercising knee flexors and thighs
- Abduction: Exercising hip and thigh muscles
- Adduction: Exercising pelvis and hips muscles
- Total Hip: Exercising hip, knee and thigh muscles
- Pull Down: Exercising arm and shoulder girdle muscles
- Shoulder Press: Exercising shoulder and arm muscles
- Chest Press: Exercising chest muscles and arm extensors
- Rowing: Exercising shoulder girdle and back muscles
- Triceps Dips: Exercising elbow extensors
- Fly: Exercising chest muscles
- Abdominal Trainer: Exercising abdominal muscles
- Low Back Trainer: Exercising back and hip muscles
- Upper Body Rotator: Exercising abdominal and back muscles

6 *Contra-indications*

6.1 **Absolute contra-indications:**

- Myocardial dysfunction in an instable phase
- HNP in the (sub) acute phase (with irradiation and neurological signs)
- Severe osteoporosis
- Acute fractures
- Severe mechanical instability
- Early post-surgical joint modification

6.2 **Relative contra-indications:**

- Myocardial disfunction in a stable phase
- HNP in the sub-acute phase (with irradiation and neurological signs)
- Mild osteoporosis
- Late post-surgical joint modification
- Tumours
- Inflammations in the involved area
- Non desired articular load in the involved area (p.e.. with M. Osgood-Schlatter, bursites, etc.)
- Non desired ligamentous load in the involved area (p.e. sprain with hydrops)
- Non desired soft tissue load in the involved area (p.e. with severe atrophy)
- Cardiovascular restrictions
- Pulmonary restrictions
- Psychological restriction
- Pain of such cause that movement would not relief

7 *Operation*

7.1 **EN-Dynamic user interface (with track)**


The training equipment in the EN-Dynamic line have uniform operation. For the EN-Dynamic Track line the operating panel is showed as fig. 1 on page 76.

1. Orange ON/OFF button.

2. Display for resistance level: This display displays the set-up resistance level during set-up and training. If a programmed Chip card with more training units is used and the last Series for this unit is finished, this display will show the next sequence number or -PC-.
Unit of resistance level [N].
3. 10 [N] = 1 [kgf].
4. “-” key for decreasing the resistance level. When changing the adjustments, the displays will flash. When the changes have been completed, you will hear a peep tone. Only start with moving after this peep tone.
5. “+” key for increasing the resistance level. When changing the adjustments, the displays will flash. When the changes have been completed, you will hear a peep tone. Only start with moving after this peep tone.
6. Next key: Pressing this key will write the data of e.g. an 1RM test to the Chip card. This function only will operate when a 1RM Chip card is loaded. This key will display the sequence numbers left in a sequential order.
In case of more than one exercises on one unit is possible (p.e. left and right) than the “>” key is also used to “accept the selected exercise”
7. Display for time indication: This display shows the duration of the training alternated with the duration of the pause between Series. Duration of the training is only displayed when this parameter is selected instead of number of Repetitions.
8. Unit of Time, (seconds)
9. Display for the Number of Series and Number of Repetitions. It counts up if time is the target. The repetition counter counts down, if “Number of Repetitions” is entered during set-up, if the counter reaches “0” an audible signal sounds. If the session is continued after reaching “0” then the counter will switch over to counting upwards.
10. “#Rep.” Unit of Repetitions and “S” number of Series.
11. “-” key will decrease the number of repetitions during set-up phase.
12. “+” key will increase the number of repetitions during the set-up phase.
13. “R” key: Pulsing this key will end the current Series and reset the Repetition counter to “0”. If more Series are programmed, the program of the training unit will enter the Pause phase before the next Series. If the Reset key is pressed during the last Series for the this unit, the display for the resistance level will show the sequence number of the next training unit on the Chip card. If pressed the “R” key when all entered sessions have finished, the resistance level display will show “-PC-” indicating that one have to return with the Chip card to the personal computer.
14. Chip card LED: If the LED lights, this indicates that the Chip card is in use performing a read or write operation to/from the card. Wait with removal of the chipcard until this LED light extinguishes.

7.2 General hardware operation

Many parts of the EN-Dynamic equipment can be changed in position. They are fixed by means of a fixation button. Be sure that the fixation button is well fixed. This means that the spring-button is not extended from the knob (see picture).

7.3 Seated Leg Press

Preparation:

- be seated with your lower back fully supported by the back rest
- adapt the foot support in such a way that your knees are approximately 90° flexed in the starting position
- place your feet on the foot support and grasp the handles

Execution:

- contract your abdominal and low back muscles before you start with the exercise
- extend your legs in a slow and controlled manner so the leg extension becomes a smooth and continuous movement

- keep the extended position for two seconds and do not overextend your knees
- flex your legs at a controlled and smooth pace

Attention marks:

- perform a complete extension and flexion
- keep the pace controlled and don't move fast to prevent any jerking movements
- keep your lower back well supported by the back rest to prevent high spinal loading

7.4 Squat

Preparation:

- Adjust the backrest to the appropriate position (sitting – lying or in-between) depending on the desired and realizable ankle- hip- and/or knee-angle
- Set the starting position of the chair (horizontal displacement)
- be seated with your lower back fully supported by the back rest and/or with your shoulders against the shoulder pads
- adapt the foot support in such a way that your knees are approximately 90° flexed in the starting position or in an alternative desired position for ankle exercising
- place your feet on the foot support and grasp the handles

Execution:

- contract your abdominal and low back muscles before you start with the exercise
- extend your legs in a slow and controlled manner so the leg extension becomes a smooth and continuous movement
- keep the extended position for two seconds and do not overextend your knees
- flex your legs at a controlled and smooth pace
- in case of an explosive movement be assured that there is ample space at the headside of the squat
- in case of an explosive movement it is recommended to start slowly in order to gain proper proprioception on launch and descent. Upon proper motor ability one can increase to full extension with a large instant of floating.

Attention marks:

- perform a complete extension and flexion
- keep the pace controlled
- keep your lower back well supported by the back rest to prevent high spinal loading

7.5 Leg Extension

Preparation:

- adapt the back angle for optimal lumbar support
- position the leg roll just above your ankle on top of your shin bone
- Set the desired starting position.
- grasp the handles

Execution:

- contract your abdominal and lower back muscles before you start with the exercise
- extend your legs in a slow and controlled manner
- keep this position for two seconds and return to the starting position with a slow and fluent movement

Attention marks:

- perform a complete extension and flexion
- do not 'lock' your knees
- control the pace and do not move too fast
- keep your shoulders suppressed and your back fully supported by the back rest

7.6 Seated Leg Curl

Preparation:

- adapt the back angle for optimal lumbar support
- position the roll below your calves just above your ankle
- Set the desired starting position.
- grasp the handles

Execution:

- contract your abdominal and lower back muscles before you start with the exercise
- bend your knees slowly at a constant pace
- keep the flexed position for two seconds and return in a constant motion, returning to the starting position

Attention marks:

- perform a complete flexion and extension
- do not 'lock' your knees
- control the pace and do not move too fast
- keep your lower back well supported by the back rest and keep your shoulders pressed down

7.7 Abduction

Preparation:

- be seated with your back fully supported by the back rest
- select the foot support so your upper legs are almost horizontal
- grasp the handles

Execution:

- contract your abdominal and lower back muscles before starting the exercise
- push your upper legs in a fluent movement outwards
- keep this position for two seconds and return to the starting position in a fluent movement

Attention marks:

- control the pace and do not move too fast
- push with your upper legs, not with your feet
- push symmetrically with both legs
- on this device one can also train unilateral (independent per leg)

7.8 Adduction

Preparation:

- be seated with your back fully supported by the back rest
- adjust the foot support so your upper legs are almost horizontal
- set the range of movement with the 2 handles right below your seat
- grasp the handles

Execution:

- contract your abdominal and lower back muscles before starting the exercise
- pull knees and thighs together in a smooth and fluent movement
- keep this position for two seconds and return to starting position with a slow and fluent movement

Attention marks:

- adjust the range of motion to a comfortable one
- control the pace and do not move too fast
- pull symmetrically with both legs
- on this device one can also train unilateral (independent per leg)

7.9 Total Hip

Preparation:

- adjust the height of the plateau so you can hold the black bar with slightly bent elbows
- position the roll at knee height
- position the roll at the front, back, inner or outer side of your knee to bring variation in the exercise
- select the desired starting angle by pulling the black button on top of the lever arm

Execution:

- contract abdominal and lower back muscles before starting the exercise
- make a fluent movement
- keep this final position for two seconds and return to the starting position with a fluent movement
- perform the exercise with both right and left legs

Attention marks:

- control the pace and do not move too fast
- keep the neck muscles relaxed and the shoulders pressed down

7.10 Pull Down

Preparation:

- select a good and comfortable seat height: knees approximately 90° flexed
- be seated with your lower back fully supported by the back rest
- place your feet on the black support

Execution:

- contract your abdominal and low back muscles before you start to exercise
- grasp the handles and pull them in a fluent pace downwards
- keep the lowest position for two seconds and move back to the starting position
- change the position of your hands focusing on different muscles (this will give variation)

Attention marks:

- perform a complete extension and flexion
- control the pace and do not move too fast
- keep your lower back well supported by the back rest and keep your neck in a relaxed position

7.11 Shoulder Press

Preparation:

- select a good and comfortable seat height: knees at approximately 90° flexed
- be seated with your lower back fully supported by the back rest
- place your feet on the foot support

Execution:

- contract your abdominal and low back muscles before you start the exercise
- take the grips and extend your arm in a gender and fluent pace
- keep this position for two seconds and return to the starting position in a slow and fluent manner

Attention marks:

- perform a complete extension and flexion
- do not 'lock' your elbows
- control the pace and do not move too fast
- keep your lower back well supported by the back rest and your neck muscles relaxed

7.12 Chest Press:

Preparation:

- select a good and comfortable seat height: knees circa 90° flexed
- be seated with your lower back fully supported by the back rest
- place your feet on the black support

Execution:

- start in a seated position and have your abdominal and back muscles contracted before you start the exercise
- grasp the handles and extend your arms away from your chest fully at a slow and fluent pace
- keep the extended position for two seconds and return to the starting position with a slow and fluent movement
- by changing the position of your hands you can concentrate attention on different parts of the muscle group

Attention marks:

- do not move your trunk or head during the exercise
- keep your shoulders pressed down in a relaxed manner
- do not 'lock' your elbows when extending your arms

7.13 Rowing

Preparation:

- adjust the seat height so your chest can comfortably rest against the chest support
- adjust the starting position by adjusting the thoracic support in such a way that your elbows are ightly flexed at the starting position
- place your feet on the support

Execution:

- contract your abdominal and lower back muscles before you start the exercise
- pull the arm of the unit to the chest in a fluent movement
- keep this position for two seconds and return to the starting position with a slow and fluent movement
- by changing the position of your hands you alter the exercise

Attention marks:

- perform a complete flexion and extension
- keep your trunk stiff by contracting abdominal and back muscles to prevent lumbar overload
- do not 'lock' your elbows and keep your shoulders and neck muscles relaxed
- control the pace and do not move too fast

7.14 Triceps Dips

Preparation:

- select a comfortable seat height: knees about 90° flexed
- be seated with your back fully supported by the back rest
- place your feet on the foot support

Execution:

- contract your abdominal and low back muscles before you start with the exercise
- take the handle in both hands (at shoulder width) and extend your elbows smoothly so that your forearms are vertical
- keep this position for two seconds and return to the starting position with a slow and fluent movement
- changing the position of your hands brings variation in the exercise

Attention marks:

- perform a complete extension and flexion
- do not overextend elbows below the horizontal and keep your shoulders pressed down
- control the pace and do not move too fast
- keep your back well supported by the back rest and your neck muscles relaxed

7.15 Fly

Preparation:

- select a good and comfortable seat height: both knees at approximately 90° flexed
- be seated with your back fully supported by the back rest
- place your feet on the foot support
- select hand grips so your upper arms are horizontal
- adjust the position of the elbow support

Execution:

- contract your abdominal and lower back muscles before you start with the exercise
- push your upper arms in a fluent pace inwards
- make a complete movement so the pads touch each other
- keep this position for two seconds and return to the starting position with a slow and fluent movement

Attention marks:

- do not only push with your hands, but with your under arms as well
- control the pace and do not move too fast
- push symmetrically
- keep the neck muscles relaxed and the shoulders pressed down
- It is also possible to exercise unilateral

7.16 Abdominal Trainer

Preparation:

- select the seat height so that your pelvis is at the same level as the unit's centre of rotation
- place the chest bar to the chest bone at a comfortable height
- have the knees bent at approximately 110° by changing the foot support
- place your feet on the foot support and below the pollster rolls

Execution:

- keep your hands crossed in front of your chest or hold the chest bar
- contract your abdominal and low back muscles before starting the exercise

- bend fluently at a slow pace forward
- keep this position for two seconds and return to the starting position with a slow and fluent movement
- for variation you can slip your feet behind the roll under the seat

Attention marks:

- keep your trunk stiff by contracting abdominal and back muscles to prevent lumbar overload

7.17 Low Back Trainer

Preparation:

- select the seat height so that your pelvis is at the same level as the unit's centre of rotation
- position the back support below your shoulder blades
- have your knees bent at approximately 110° by changing the foot support
- place your feet on the foot support and below the pollster rolls

Execution:

- keep your arms crossed in front of your chest
- contract your abdominal and low back muscles before starting the exercise
- extend your back slowly at a fluent pace
- keep the extended position for two seconds and return to the starting position with a slow and fluent movement

Attention marks:

- do not move too fast and do not over-extend your back
- keep your trunk stiff by contracting abdominal and back muscles to prevent lumbar overload
- have the back support positioned on a comfortable height to perform the exercise with the right force

7.18 Upper Body Rotator

Preparation:

- set the starting position of the arm

Execution:

- turn clock-wise: place the left-arm behind the left-roll and the right arm in front of the right roll
- turn counter clock-wise: place the right-arm behind the right roll and the left arm in front of the left roll
- contract your abdominal and low back muscles before starting the exercise
- turn with the upper-body back to the neutral starting position
- keep the extended position for two seconds and return to the starting position with a slow and fluent movement

Attention marks:

- rotate the head in the same direction as the body
- keep your trunk stiff by contracting abdominal and back muscles to prevent lumbar overload

8 Cleaning and maintenance

For correct operation keep the bars of the mechanism clean and free of grease.

The EN-Dynamic may be cleaned with a damp cloth and a domestic liquid cleaner. Use lukewarm water and a non-abrasive liquid household cleaner (no abrasive, no alcohol content solution).

8.1 Daily maintenance

- Check all moving parts with regard to wear or visible deviations such as damage and cracks
- Check the belts, if visible, with regard to wear
- Check the pneumatic pipes (auditive) with report to possible leakage

8.2 Weekly maintenance

- Remove any condensation moisture from the compressor (consult the operation manual of your compressor for this).

8.3 Periodical maintenance

Periodical maintenance depends on the intensity of use. With an average use (8 hours a day, 5 days a week) we advise you to perform the following maintenance periodically (every month).

- Clean the piston rods
- If necessary, lubricate the points of rotation with light grade oil.

8.4 Annual maintenance

- Let an authorised Enraf-Nonius service employee check your equipment every year

In the event of possible failures or deviations observed by you during one of these check-ups, you must take care that the equipment cannot be used anymore until these failures or deviations have been completely solved.

9 Finally

Your EN-Dynamic contains materials which can be recycled and/or are noxious for the environment. At the end of the duration of life, specialised companies can take apart the EN-Dynamic and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment. Please ensure that you are well informed of the local rules and regulations regarding to the removal of equipment and accessories.

We hope that you will use the EN-Dynamic for a long time with much satisfaction and we take the liberty of drawing your attention to the many other high quality products in our programme. If you should have any questions or wish to receive our catalogue, please contact our representative or visit us at www.enraf-nonius.com.

10 Technical data

1455913 EN-Dynamic Seated Leg Press (external low voltage power adapter)

1455933 EN-Dynamic Seated Leg Press (internal power supply with mains cable)

- Range: 0 - 200 x 10 N (kgf)
- Size: 220 x 90 x 140 cm
- Weight: 105 kg
- Transport information: 183 x 98 x 161 cm/135 kg

1455915 EN-Dynamic Squat (external low voltage power adapter)

1455935 EN-Dynamic Squat (internal power supply with mains cable)

- Range: 0 - 200 x 10 N (kgf)
- Size: 250 x 80 x 150 cm
- Weight: 265 kg
- Transport information: 280 x 100 x 160 cm/280 kg

1455905 EN-Dynamic Leg Extension (external low voltage power adapter)

1455925 EN-Dynamic Leg Extension (internal power supply with mains cable)

- Range: 0 - 100 x 10 N (kgf)
- Size: 140 x 90 x 140 cm
- Weight: 105 kg
- Transport information: 148 x 88 x 161 cm/135 kg

1455912 EN-Dynamic Seated Leg Curl (external low voltage power adapter)

1455932 EN-Dynamic Seated Leg Curl (internal power supply with mains cable)

- Range: 0 - 70 x 10 N (kgf)
- Size: 140 x 90 x 140 cm
- Weight: 105 kg
- Transport information: 128 x 103 x 161 cm/135 kg

1455909 EN-Dynamic Abduction (external low voltage power adapter)

1455929 EN-Dynamic Abduction (internal power supply with mains cable)

- Range: 0 - 60 x 10 N (kgf)
- Size: 140 x 140 x 140 cm
- Weight: 85 kg

- Transport information: 148 x 88 x 161 cm/115 kg

1455908 EN-Dynamic Adduction (external low voltage power adapter)

1455928 EN-Dynamic Adduction (internal power supply with mains cable)

- Range: 0 - 60 x 10 N (kgf)
- Size: 140 x 140 x 140 cm
- Weight: 85 kg
- Transport information: 148 x 88 x 161 cm/115 kg

1455911 EN-Dynamic Total Hip (external low voltage power adapter)

1455931 EN-Dynamic Total Hip (internal power supply with mains cable)

- Range: 0 - 36 x 10 N (kgf)
- Size: 110 x 115 x 170 cm
- Weight: 120 kg
- Transport information: 128 x 103 x 161 cm/150 kg

1455903 EN-Dynamic Pull Down (external low voltage power adapter)

1455923 EN-Dynamic Pull Down (internal power supply with mains cable)

- Range: 0 - 80 x 10 N (kgf)
- Size: 140 x 120 x 180 cm
- Weight: 105 kg
- Transport information: 158 x 128 x 161 cm/135 kg

1455907 EN-Dynamic Shoulder Press (external low voltage power adapter)

1455927 EN-Dynamic Shoulder Press (internal power supply with mains cable)

- Range: 0 - 60 x 10 N (kgf)
- Size: 140 x 120 x 180 cm
- Weight: 135 kg
- Transport information: 188 x 128 x 161 cm/165 kg

1455901 EN-Dynamic Chest Press (external low voltage power adapter)

1455921 EN-Dynamic Chest Press (internal power supply with mains cable)

- Range: 0 - 120 x 10 N (kgf)
- Size: 135 x 105 x 140 cm
- Weight: 100 kg
- Transport information: 158 x 128 x 161 cm/130 kg

1455906 EN-Dynamic Rowing (external low voltage power adapter)

1455926 EN-Dynamic Rowing (internal power supply with mains cable)

- Range: 0 - 100 x 10 N (kgf)
- Size: 120 x 90 x 115 cm
- Weight: 100 kg
- Transport information: 128 x 103 x 136 cm/130 kg

1455914 EN-Dynamic Triceps Dips (external low voltage power adapter)

1455934 EN-Dynamic Triceps Dips (internal power supply with mains cable)

- Range: 0 - 100 x 10 N (kgf)
- Size: 130 x 72 x 140 cm
- Weight: 96 kg
- Transport information: 148 x 88 x 161 cm/126 kg

1455910 EN-Dynamic Fly (external low voltage power adapter)

1455930 EN-Dynamic Fly (internal power supply with mains cable)

- Range: 0 - 65 x 10 N (kgf)
- Size: 130 x 120 x 140 cm
- Weight: 100 kg
- Transport information: 158 x 128 x 161 cm/130 kg

1455904 EN-Dynamic Abdominal Trainer (external low voltage power adapter)

1455924 EN-Dynamic Abdominal Trainer (internal power supply with mains cable)

- Range: 0 - 80 x 10 N (kgf)
- Size: 160 x 90 x 140 cm
- Weight: 120 kg
- Transport information: 183 x 98 x 161 cm/150 kg

GB

27

1455902 EN-Dynamic Low Back Trainer (external low voltage power adapter)
1455922 EN-Dynamic Low Back Trainer (internal power supply with mains cable)

- Range: 0 - 80 x 10 N (kgf)
- Size: 160 x 90 x 130 cm
- Weight: 120 kg
- Transport information: 183 x 98 x 161 cm/140 kg

1455936 EN-Dynamic Rotator (external low voltage power adapter)

- Range: 0 - 60 x 10 N (kgf)
- Size: 135 x 95 x 108 cm
- Weight: 85 kg
- Transport information: 153 x 105 x 128 cm/115 kg

In case of a connection of an external adapter, other than the prescribed charger can adversely affect the safety of the user and the functioning of the unit, and is therefore not permitted. Allowed are:

- 3444488 Adaptor 230V EUR
- 3444489 Adapter 230/12V - BSI (GB) for EN-Dynamic and EN-Train

In case of malfunction the actual resistance may deviate 20% from displayed resistance.

Environment conditions for transport and storage

Environment temperature : -10° till +50° C
Relative humidity : 10 till 95 %
Atmospheric pressure : 500 till 1060 hPa

Environments conditions for normal use:

Environment temperature : 10° till 40° C
Relative humidity : 10 till 90 %
Atmospheric pressure : 500 till 1060 hPa

10.1 Power supply

External low voltage power adapter

Input : 230 V / 50-60 Hz
Output : 12V – 10W

Internal power supply with mains cable

Power supply : 90-240V / 50-60Hz / 15VA

10.2 Classification

This equipment complies with all requirements of the Medical Devices Directive (93/42/EEC).

Medical device classification : Klasse IIa

Safety class (IEC 60601-1) (IEC 60601-1)

External low voltage power adapter : II (double insulated)

Internal power supply with mains cable: Klasse I B

10.3 Installation requirements EN-Dynamic

The installation requirements for EN-Dynamic are:

- accessibility of the room, i.e.
- doors of at least 85 cm width
- preferably no stairs (units weigh up to 265 kg !!)
- when elevators are being used, check the doorway, the weight capability and the sizes of the elevator (a Squat is 250 cm long!)

For the EN-Dynamic installation one needs

- enough wall outlets, one for every EN-Dynamic unit and one for the compressor
- a skirtboard where the air tubes can be put in (p.e. tehait boards. The air tubes have an 8 mm diameter)

The compressor is whispering silent, thus it can be installed in the same room as the EN-Dynamic units. It however can also be installed in a separate room.

If the institute has it's own air supply system, assure there is 8,5 - 9,5 bar pressure. Less pressure (up to 6 bar) will result in longer setting time of the units, and the maximal performance (maximal resistance setting) cannot be achieved. The compressor provided by Enraf-Nonius is a Jun-Air Model "6-15-Standard").

In case of a locally obtained compressor, please consult the Jun-air website for your local dealer (<http://www.jun-air.dk>). The specifications for the compressor are:

- 15 litre pressure barrel (air receiver)
- Cut-in pressure: 8,5 bar
- Cut-out pressure: 9,5 bar
- included a pressure gauge and an outlet cock
- 0,5mu moisture filter and 0,01mu oil filter
- connector 1/4", 230V 50Hz

10.4 Installation requirements EN-Track

Installation requirements are equal to the EN-Dynamic requirements, plus a computer and cardreader to run the EN-Track software and to write and read the chipcards.

EN-Track PC configuration:

- Operating system: Windows XP
- Hardware requirement: cd-rom drive, (S)VGA monitor, USB interface
- EN-Track hosting PC's should be cable connected (no WLAN)

Included in the Track installation kit are 10 chipcards.

They can be reloaded for every patient, and thus one needs only as much chipcards as one has exercise systems. In case one wants to provide every user with his own private card, then additional cards can be ordered (in packages of 10).

10.5 Guidance and manufacturer 's declaration –electromagnetic immunity

The EN-Dynamic is intended for use in the electromagnetic environment specified below. The customer or the user of the EN-Dynamic should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV kV contact ± 8 kV air	± 6 kV kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 1 kV for power supply lines not applicable	Mains power quality should be that of a typical computer room.
Surge IEC 61000-4-5	± 1 kV line(s)to line(s) ± 2 kV line(s)to earth	± 1 kV line(s)to line(s) ± 2 kV line(s)to earth	Mains power quality should be that of a typical computer room.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	< 5% U_T for 0,5 cycle < 5% U_T for 1 cycle 70% U_T for 25 cycles < 5% U_T for 5 sec.	Mains power quality should be that of a typical computer room. If the user of the EN-Dynamic requires continued operation during power mains interruptions, it is recommended that the EN-Dynamic be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3A / m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical computer room.

NOTE U_T is the a.c.mains voltage prior to application of the test level.

10.6 Guidance and manufacturer's declaration –electromagnetic immunity

The EN-Dynamic is intended for use in the electromagnetic environment specified below. The customer or the user of the EN-Dynamic should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the EN-Dynamic, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [3,5/3]\sqrt{P}$ $d = [3,5/3]\sqrt{P}$ 80 MHz to 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/ m 80 MHz to 2,5 GHz	3 V/ m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the EN-Dynamic is used exceeds the applicable RF compliance level above, the EN-Dynamic should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the EN-Dynamic.

b

Over the frequency range 15 kHz to 8 MHz, field strengths should be less than 3 V / m.

10.7 Recommended separation distances between portable and mobile RF communications equipment and the EN-Dynamic

The EN-Dynamic is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the EN-Dynamic can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the EN-Dynamic as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = [3,5/V_1]\sqrt{P}$	80 MHz to 800 MHz $d = [3,5/E_1]\sqrt{P}$	800 MHz to 2,5 GHz $d = [7/E_1]\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,69	3,69	7,38
100	11,67	11,67	23,33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

11 Ordering data

For the ordering data of the EN-Dynamic, standard accessories and additional accessories we refer to the Catalogue for Physiotherapy and Rehabilitation or to www.enraf-nonius.com.

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